A1CNow⁺ Controls™ Levels 1+2

Svenska

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Analytes / Instruments	Units	Level 1 / Lot	Level 2 / Lot	SI Units¹	Level 1 / Lot	Level 2 / Lot
	mean	expected range	mean	expected range	mean	expected range
Hemoglobin A1c						
A1CNow ⁺ System	%			mmol/mol		

Footnotes for values | Fußnoten für werte | Apostilles pour des valeurs | Note a piè di pagina per i valori | Notas al pie de la página para los valores
Voetnoten voor waarden | Notas de rodapé para valores | Przepisy dotyczące wartości | Fotnoter till värden

English	Français	Español	Portuguese	Svenska
1 SI International System of Units	1 Système international d’unités, SI	1 SI, Sistema Internacional de unidades	1 Sistema Internacional de Unidades (SI)	1 SI Internationella måttenhetssystemet
Deutsch	Italiano	Nederlands	Polskie	
1 Internationales Einheitensystem	1 Sistema internazionale di unità di misura	1 Internationale Stelsel van Eenheden	1 Międzynarodowy system jednostek	

Polymer Technology Systems, Inc.
4600 Anson Boulevard
Whitestown, IN 46075 USA
+1-317-870-5610

PTS Diagnostics 4600 Anson Boulevard, Whitestown, IN 46075 USA phone: (317) 870-5610 fax (317) 870-5608 www.ptsdiagnostics.com

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English

Intended Use

The A1CNow⁺ Controls™ are intended for the quality control of PTS Diagnostics' A1CNow HbA1c tests.

Product Description

The A1CNow⁺ Controls are supplied in two levels as a frozen, ready-to-use material that requires no reconstitution. The control consists of human whole blood modified to target levels with reagent chemicals to produce the desired reaction when tested by the A1CNow system. Stabilizers and preservatives have been added to maintain shelf life and inhibit microbial growth.

Warnings and Precautions

POTENTIAL BIOHAZARDOUS MATERIAL. All blood donor units used in this product have been tested and found non-reactive for Hepatitis B Surface Antigen, Hepatitis C and HIV 1 & 2 antibody when tested by FDA accepted methods.

No known test method can assure that a product derived from human blood does not contain Hepatitis or HIV viruses. Handle QC material as you would a patient sample. QC materials should be used and disposed of in accordance with regulatory and accreditation requirements.

Storage and Stability

The A1CNow⁺ Controls should be stored frozen at -30 to -10°C (-22 to 14°F). The control may be received thawed and should be re-frozen upon arrival. When stored frozen at -30 to -10°C (-22 to 14°F), the control remains stable until the expiration date printed on the label. Once the control is thawed and opened, it may be stored refrigerated at 2-8°C (36-46°F) for 8 months or at room temperature at 18 to 25°C (64 to 77°F) for 7 days. Discard the control if there is any evidence of microbial contamination. Discard controls in the same manner as other biological specimens according to local guidelines.

Procedure

If frozen, allow the control to thaw for at least 30 minutes at room temperature with the cap on. Mix gently by inversion a minimum of ten times to ensure the homogeneity of the contents. Avoid foaming. Remove cap and squeeze the proper amount of control into a sample cup or on a clean flat surface for assay. Treat the control as you would a patient sample in accordance with the manufacturer's requirements for the test method. Immediately recap and return the control to its desired storage condition. Maintain the same storage condition for the life of the product. Alternating the control between refrigerated temperature and room temperature storage could adversely affect the product's shelf life. If the control is received thawed, immediately place in the desired closed vial storage condition. For further information see the storage and stability section of this insert.

Expected Values

Expected values for the listed lots of controls have been established from inter-laboratory data using the A1CNow system. Individual laboratory means should fall within the ranges listed. These values should be used as a guide in evaluating the performance of test methods. Each laboratory should establish its own ranges and specifications for the methods used to measure HbA1c.

Limitations

The expected mean and ranges were established using the A1CNow system available at the time of assay. Any future changes made to the test method may give different values from those previously recovered. Use of methods other than the A1CNow system may give different values from the ones indicated. Limitations of the test method are included in the package insert for the A1CNow system. Technical updates can be found on our website.

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