A1CNow⁺ Controls™ Levels 1+2

LOT 61041
61042

REF 737

2022-11-30

CE IVD

-30°C

CONTROL

European Conformity
CE-Konformitätskennzeichnung
Conformité aux normes
européennes
Conformità europea
Conformidade europea
Europäische Konformität
Conformidade Europeia
Згідність з вимогами ЄС
EU-verenstemming

Lot Number
Designation du lot
Numero di lotto
Denominación de lote
Partijnummer
Nº do lote
Numer partii
Partinummer

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English

Intended Use

The A1CNow⁺ Controls™ are intended for the quality control of PTS Diagnostics' A1CNow HbA1c tests.

Product Description

The A1CNow⁺ Controls are supplied in two levels as a frozen, ready-to-use material that requires no reconstitution. The control consists of human whole blood modified to target levels with reagent chemicals to produce the desired reaction when tested by the A1CNow system. Stabilizers and preservatives have been added to maintain shelf life and inhibit microbial growth.

Warnings and Precautions

POTENTIAL BIOHAZARDOUS MATERIAL. All blood donor units used in this product have been tested and found non-reactive for Hepatitis B Surface Antigen, Hepatitis C and HIV 1 & 2 antibody when tested by FDA accepted methods.

No known test method can assure that a product derived from human blood does not contain Hepatitis or HIV viruses. Handle QC material as you would a patient sample. QC materials should be used and disposed of in accordance with regulatory and accreditation requirements.

Storage and Stability

The A1CNow⁺ Controls should be stored frozen at -30 to -10°C (-22 to 14°F). The control may be received thawed and should be re-frozen upon arrival. When stored frozen at -30 to -10°C (-22 to 14°F), the control remains stable until the expiration date printed on the label. Once the control is thawed and opened, it may be stored refrigerated at 2-8°C (36-46°F) for 8 months or at room temperature at 18 to 25°C (64 to 77°F) for 7 days. Discard the control if there is any evidence of microbial contamination. Discard controls in the same manner as other biological specimens according to local guidelines.

Procedure

If frozen, allow the control to thaw for at least 30 minutes at room temperature with the cap on. Mix gently by inversion a minimum of ten times to ensure the homogeneity of the contents. Avoid foaming. Remove cap and squeeze the proper amount of control into a sample cup or on a clean flat surface for assay. Treat the control as you would a patient sample in accordance with the manufacturer's requirements for the test method. Immediately recap and return the control to its desired storage condition. Maintain the same storage condition for the life of the product. Alternating the control between refrigerated temperature and room temperature storage could adversely affect the product's shelf life. If the control is received thawed, immediately place in the desired closed vial storage condition. For further information see the storage and stability section of this insert.

Expected Values

Expected values for the listed lots of controls have been established from inter-laboratory data using the A1CNow system. Individual laboratory means should fall within the ranges listed. These values should be used as a guide in evaluating the performance of test methods. Each laboratory should establish its own ranges and specifications for the methods used to measure HbA1c.

Limitations

The expected mean and ranges were established using the A1CNow system available at the time of assay. Any future changes made to the test method may give different values from those previously recovered. Use of methods other than the A1CNow system may give different values from the ones indicated. Limitations of the test method are included in the package insert for the A1CNow system. Technical updates can be found on our website.

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Polymer Technology Systems, Inc.
MDSS GmbH