

# ChekMate Strips Range Card

Please enter ranges listed below on the Results Sheet included with your ChekMate kit.

Batch #:

0860-31017-01

Expiration date:

2019 05 06

	C1 RED	C2 RED	C3 RED	C1 GREEN
Level 1	21.50 – 28.00	21.50 – 28.00	20.00 – 27.00	15.00 – 23.00
Level 2	37.00 – 45.00	37.50 – 45.00	37.75 – 43.75	42.75 – 51.25

ChekMate results are reported as % reflectance and should not be compared to values obtained with PTS Panels® test strips.

PTS Diagnostics Customer Service: 877-870-5610 | Fax: 317-870-5608 | Outside of US: +1-317-870-5610

Email: [customerservice@ptsdiagnostics.com](mailto:customerservice@ptsdiagnostics.com)

PS-000882 Rev. 2