

Does Precise Sample Size Affect Your Outcome?

Bulb-squeeze dispensing technology ensures consistent performance

Plastic casing decreases breakage

Glass tube provides consistent capillary action

No assembly required



Ideal for all blood collection needs up to 50 µL

Speeds Blood Collection and Ensures Savings

PTS Collect™ capillary tubes provide significant savings over traditional transfer tube technology—especially for applications where exact volumes are required. Using a single glass tube with plastic casing, PTS Collect capillary tubes provide consistent self-wicking technology to preset volumes.

Increased Efficiency in Workflow

- Ease-of-use with no assembly
- One-hand blood collection
- Glass tube ensures easy blood collection
- Bulb-squeeze dispensing technology ensures consistent performance
- Only one item to inventory (no plunger)

Sample Volume	15 µL	20 µL	30 µL	40 µL	50 µL
Quantity	25	25	25	16	25
Ref. #	[2863]	[2864]	[2865]	[2866]	[2134]



PTS Diagnostics
 7736 Zionsville Road, Indianapolis, IN 46268 USA
 +1 (317) 870-5610