


# PTS Collect™ Capillary Tubes Consistently Collect Precise Blood Samples<sup>1</sup>



PTS Collect Capillary Tube 40 µL Sample Volume  
Results for all available sizes fall between low and high tolerances



## Advanced Technology for Point-of-Care Testing

PTS Collect capillary tubes utilize a single glass tube with plastic casing and provide consistent self-wicking technology to preset volumes. This innovative technology is more efficient, offers significant savings over traditional blood transfer tubes, and is designed to be easier to use.

The improved design is ideal for most small volume blood collection needs.



Sample Volume	15 µL	20 µL	30 µL	40 µL	50 µL
Quantity	25	25	25	16	25
Ref. #	[2863]	[2864]	[2865]	[2866]	[2134]



To order, contact your local authorized distributor.

PTS Diagnostics  
7736 Zionsville Road, Indianapolis, IN 46268 USA  
+1 (317) 870-5610

1. PTS Diagnostics Capillary Tubes First Article Inspection Summary, January 2015. Data on file.