

pts panels[®]

HDL Cholesterol Controls **Range Card**

Lot #: HC27

Expiration Date: 2020-04-30

Single Analyte Test Strips

HDL Cholesterol*
(*CardioChek[®] PA analyzers and CardioChek[®] Plus analyzers version 1.08 or below*)

Multi-Analyte Test Strips

HDL Cholesterol*
(*CardioChek[®] PA analyzers and CardioChek[®] Plus analyzers version 1.08 or below*)

Level 1

mg/dL

mmol/L

≤25-45

≤0.65-1.17

≤15-50

≤0.39-1.30

Level 2

mg/dL

mmol/L

30-70

0.78-1.81

30-70

0.78-1.81

* For CardioChek, CardioChek PA, and CardioChek Plus (version 1.08 or below) analyzers.

PTS Diagnostics Customer Service: 877-870-5610 | Fax: 317-870-5608 | Outside of US: +1-317-870-5610
Email: customerservice@ptsdiagnostics.com

PS-002516 Rev. 7