

<b>pts panels<sup>®</sup></b> Multi-Chemistry Controls      Range Card	Level 1		Level 2	
	mg/dL	mmoL/L	mg/dL	mmoL/L
<b>Lot #:</b> <b>MC27</b> <b>Expiration Date:</b> <b>2019 02 28</b>				
<b>Single-Analyte Test Strips</b>				
Cholesterol	120-230	3.11-5.96	145-330	3.76-8.55
HDL Cholesterol* <i>(CardioChek<sup>®</sup> Plus analyzer version 1.09 or above)</i>	≤25-50	≤0.65-1.30	30-80	0.78-2.07
Triglycerides	95-190	1.07-2.15	150-270	1.70-3.05
Glucose/eGLU	40-125	2.22-6.94	170-380	9.44-21.09
<b>Multi-Analyte Test Strips</b>				
Cholesterol	115-230	2.98-5.96	150-330	3.89-8.55
HDL Cholesterol* <i>(CardioChek<sup>®</sup> Plus analyzer version 1.09 or above)</i>	≤15-55	≤0.39-1.42	30-90	0.78-2.33
Triglycerides	70-190	0.79-2.15	140-320	1.58-3.62
Glucose	20-90	1.11-5.00	110-340	6.11-18.87
* For CardioChek, CardioChek PA, and CardioChek Plus (version 1.08 or below) analyzers, HDL Cholesterol requires a separate HDL Cholesterol Control (Cat. No. 722).				
PTS Diagnostics Customer Service: 877-870-5610   Fax: 317-870-5608   Outside of US: +1-317-870-5610 Email: customerservice@ptsdiagnostics.com				
PS-002519 Rev. 11				