

pts panels[®] Multi-Chemistry Controls Range Card	Level 1		Level 2	
	mg/dL	mmoL/L	mg/dL	mmoL/L
Lot #: MC29 Expiration Date: 2019 07 31				
Single-Analyte Test Strips				
Cholesterol	120-230	3.11-5.96	150-330	3.89-8.55
HDL Cholesterol* <i>(CardioChek[®] Plus analyzer version 1.09 or above)</i>	≤25-50	≤0.65-1.30	30-80	0.78-2.07
Triglycerides	100-200	1.13-2.26	160-280	1.81-3.16
Glucose/eGLU	40-130	2.22-7.22	170-380	9.44-21.09
Multi-Analyte Test Strips				
Cholesterol	110-230	2.85-5.96	140-330	3.63-8.55
HDL Cholesterol* <i>(CardioChek[®] Plus analyzer version 1.09 or above)</i>	≤15-60	≤0.39-1.55	30-90	0.78-2.33
Triglycerides	80-200	0.90-2.26	150-330	1.70-3.73
Glucose	20-100	1.11-5.55	120-360	6.66-19.98
* For CardioChek, CardioChek PA, and CardioChek Plus (version 1.08 or below) analyzers, HDL Cholesterol requires a separate HDL Cholesterol Control (Cat. No. 722).				
PTS Diagnostics Customer Service: 877-870-5610 Fax: 317-870-5608 Outside of US: +1-317-870-5610 Email: customerservice@ptsdiagnostics.com				
PS-002519 Rev. 11				