

<b>pts panels<sup>®</sup></b> Multi-Chemistry Controls      Range Card		Level 1		Level 2	
		mg/dL	mmoL/L	mg/dL	mmoL/L
Lot #:	<b>MC31</b>				
Expiration Date:	<b>2020 03 31</b>				
<b>Single-Analyte Test Strips</b>					
Cholesterol		120-230	3.11-5.96	150-330	3.89-8.55
HDL Cholesterol*	<i>(CardioChek<sup>®</sup> Plus analyzer version 1.09 or above)</i>	≤25-50	≤0.65-1.30	30-80	0.78-2.07
Triglycerides		100-200	1.13-2.26	180-300	2.03-3.39
Glucose/eGLU		40-130	2.22-7.22	160-370	8.88-20.54
<b>Multi-Analyte Test Strips</b>					
Cholesterol		110-230	2.85-5.96	140-330	3.63-8.55
HDL Cholesterol*	<i>(CardioChek<sup>®</sup> Plus analyzer version 1.09 or above)</i>	≤15-60	≤0.39-1.55	30-90	0.78-2.33
Triglycerides		90-200	1.02-2.26	150-320	1.70-3.62
Glucose		20-100	1.11-5.55	125-350	6.94-19.43
* For CardioChek, CardioChek PA, and CardioChek Plus (version 1.08 or below) analyzers, HDL Cholesterol requires a separate HDL Cholesterol Control (Cat. No. 722).					
PTS Diagnostics Customer Service: 877-870-5610   Fax: 317-870-5608   Outside of US: +1-317-870-5610 Email: customerservice@ptsdiagnostics.com					
PS-002519 Rev. 11					