

Customer:	State:						
Event ID:							
PTS Panels® Test Strips							
Test Strip/MEMo Chip®	Lot #:	Exp:					
Multi-Chemistry Control	Lot #:	Exp:					
HDL Cholesterol Control	Lot #:	Exp:					

Result Range	Level 1	Level 2				
Total Cholesterol						
Triglycerides						
GLU/eGLU™						
HDL Cholesterol*						

Important:

- 1. Use the gray CardioChek® Check Strips to conduct an Optics Test prior to running controls.
- 2. Do not use PTS Panels® Control Solution beyond the expiration date on the vial. Open product is stable for 10 months at room temperature.

Serial	#:				,	•	Optics	Check	□ P	ass	□ Fail						
		Level 1					Level 2				Repeat: ■ Level 1 ■ Level 2						
Date	Initial	CHOL	TRG	GLU	HDL	Pass/Fail	CHOL	TRG	GLU	HDL	Pass/Fail	CHOL	TRG	GLU	HDL	Pass/Fail	
						□P/□F					□P/□F					□P/□F	
Signature: Notes:																	
Serial #: Optics Check: Pass Fail																	
		Level 1			Level 2				Repeat: ■ Level 1 ■ Level				evel 2				
Date	Initial	CHOL	TRG	GLU	HDL	Pass/Fail	CHOL	TRG	GLU	HDL	Pass/Fail	CHOL	TRG	GLU	HDL	Pass/Fail	
						□P/□F					□P/□F					□P/□F	
Signatu	ıre:						Notes:										
Serial	#:						Optics	Check	: □P	ass	□ Fail						
		Level 1							Level	2		Rep	oeat: I	Level	1 🔲 L	evel 2	
Date	Initial	CHOL	TRG	GLU	HDL	Pass/Fail	CHOL	TRG	GLU	HDL	Pass/Fail	CHOL	TRG	GLU	HDL	Pass/Fail	
						□P/□F					□P/□F					□P/□F	
Signatu	ıre:						Notes:										
Serial	#:						Optics	s Check	:	ass	□ Fail						
				Level	1				Level	2		Repeat: ■ Level 1			1 🔲 L	Level 2	
Date	Initial	CHOL	TRG	GLU	HDL	Pass/Fail	CHOL	TRG	GLU	HDL	Pass/Fail	CHOL	TRG	GLU	HDL	Pass/Fail	
						□P/□F					□P/□F					□P/□F	
Signati	Signature: Notes:																
Serial #: Optics Check: □ Pass □ Fail																	
				Level	1				Level	2		Re	peat: I	Level	1 🔲 L	evel 2	
Date	Initial	CHOL	TRG	GLU	HDL	Pass/Fail	CHOL	TRG	GLU	HDL	Pass/Fail	CHOL	TRG	GLU	HDL	Pass/Fail	
						□P/□F					□P/□F					□P/□F	
Signati	ire.						Notes:										

CardioChek, PTS Panels, eGLU and MEMo Chip are trademarks of Polymer Technology Systems, Inc. © 2015 Polymer Technology Systems, Inc. MKG 001156 Rev.2 02/15



Problem and Corrective Action Log:

^{*}HDL Cholesterol Control is done separately