

pts panels[®] Multi-Chemistry Controls Range Card	Level 1		Level 2	
	mg/dL	mmoL/L	mg/dL	mmoL/L
Lot #: MC35 Expiration Date: 2021-04-30				
Single-Analyte Test Strips				
Cholesterol	110-220	2.85-5.70	130-310	3.37-8.03
HDL Cholesterol* <i>(CardioChek[®] Plus analyzer version 1.09 or above)</i>	≤25-60	≤0.65-1.55	30-80	0.78-2.07
Triglycerides	90-200	1.02-2.26	140-290	1.58-3.28
Glucose/eGLU	40-130	2.22-7.22	150-370	8.33-20.54
Multi-Analyte Test Strips				
Cholesterol	105-225	2.72-5.83	130-320	3.37-8.29
HDL Cholesterol* <i>(CardioChek[®] Plus analyzer version 1.09 or above)</i>	≤15-60	≤0.39-1.55	30-90	0.78-2.33
Triglycerides	70-200	0.79-2.26	140-320	1.58-3.62
Glucose	20-100	1.11-5.55	110-350	6.11-19.43
* For CardioChek, CardioChek PA, and CardioChek Plus (version 1.08 or below) analyzers, HDL Cholesterol requires a separate HDL Cholesterol Control (Cat. No. 722). PTS Diagnostics Customer Service: 877-870-5610 Fax: 317-870-5608 Outside of US: +1-317-870-5610 Email: customerservice@ptsdiagnostics.com				
PS-002519 Rev. 11				

pts panels[®] Multi-Chemistry Controls Range Card	Level 1		Level 2	
	mg/dL	mmoL/L	mg/dL	mmoL/L
Lot #: MC35 Expiration Date: 2021-04-30				
Single-Analyte Test Strips				
Cholesterol	110-220	2.85-5.70	130-310	3.37-8.03
HDL Cholesterol* <i>(CardioChek[®] Plus analyzer version 1.09 or above)</i>	≤25-60	≤0.65-1.55	30-80	0.78-2.07
Triglycerides	90-200	1.02-2.26	140-290	1.58-3.28
Glucose/eGLU	40-130	2.22-7.22	150-370	8.33-20.54
Multi-Analyte Test Strips				
Cholesterol	105-225	2.72-5.83	130-320	3.37-8.29
HDL Cholesterol* <i>(CardioChek[®] Plus analyzer version 1.09 or above)</i>	≤15-60	≤0.39-1.55	30-90	0.78-2.33
Triglycerides	70-200	0.79-2.26	140-320	1.58-3.62
Glucose	20-100	1.11-5.55	110-350	6.11-19.43
* For CardioChek, CardioChek PA, and CardioChek Plus (version 1.08 or before) analyzers, HDL Cholesterol requires a separate HDL Cholesterol Control (Cat. No. 722).. PTS Diagnostics Customer Service: 877-870-5610 Fax: 317-870-5608 Outside of US: +1-317-870-5610 Email: customerservice@ptsdiagnostics.com				
PS-002519 Rev. 11				