

<b>pts panels<sup>®</sup></b> Multi-Chemistry Controls      Range Card	Level 1		Level 2	
	mg/dL	mmoL/L	mg/dL	mmoL/L
<b>Lot #: MC36</b> <b>Expiration Date: 2021-08-31</b>				
<b>Single-Analyte Test Strips</b>				
Cholesterol	110-220	2.85-5.70	140-320	3.63-8.29
HDL Cholesterol* <i>(CardioChek<sup>®</sup> Plus analyzer version 1.09 or above)</i>	≤25-60	≤0.65-1.55	30-80	0.78-2.07
Triglycerides	90-200	1.02-2.26	150-300	1.70-3.39
Glucose/eGLU	40-130	2.22-7.22	150-380	8.33-21.09
<b>Multi-Analyte Test Strips</b>				
Cholesterol	105-225	2.72-5.83	130-320	3.37-8.29
HDL Cholesterol* <i>(CardioChek<sup>®</sup> Plus analyzer version 1.09 or above)</i>	≤15-60	≤0.39-1.55	25-85	0.65-2.20
Triglycerides	70-200	0.79-2.26	140-320	1.58-3.62
Glucose	20-110	1.11-6.11	115-360	6.38-19.98
* For CardioChek, CardioChek PA, and CardioChek Plus (version 1.08 or below) analyzers, HDL Cholesterol requires a separate HDL Cholesterol Control (Cat. No. 722).				
PTS Diagnostics Customer Service: 877-870-5610   Fax: 317-870-5608   Outside of US: +1-317-870-5610 Email: customerservice@ptsdiagnostics.com				
PS-002519 Rev. 11				

<b>pts panels<sup>®</sup></b> Multi-Chemistry Controls      Range Card	Level 1		Level 2	
	mg/dL	mmoL/L	mg/dL	mmoL/L
<b>Lot #: MC36</b> <b>Expiration Date: 2021-08-31</b>				
<b>Single-Analyte Test Strips</b>				
Cholesterol	110-220	2.85-5.70	140-320	3.63-8.29
HDL Cholesterol* <i>(CardioChek<sup>®</sup> Plus analyzer version 1.09 or above)</i>	≤25-60	≤0.65-1.55	30-80	0.78-2.07
Triglycerides	90-200	1.02-2.26	150-300	1.70-3.39
Glucose/eGLU	40-130	2.22-7.22	150-380	8.33-21.09
<b>Multi-Analyte Test Strips</b>				
Cholesterol	105-225	2.72-5.83	130-320	3.37-8.29
HDL Cholesterol* <i>(CardioChek<sup>®</sup> Plus analyzer version 1.09 or above)</i>	≤15-60	≤0.39-1.55	25-85	0.65-2.20
Triglycerides	70-200	0.79-2.26	140-320	1.58-3.62
Glucose	20-110	1.11-6.11	115-360	6.38-19.98
* For CardioChek, CardioChek PA, and CardioChek Plus (version 1.08 or before) analyzers, HDL Cholesterol requires a separate HDL Cholesterol Control (Cat. No. 722).				
PTS Diagnostics Customer Service: 877-870-5610   Fax: 317-870-5608   Outside of US: +1-317-870-5610 Email: customerservice@ptsdiagnostics.com				
PS-002519 Rev. 11				