

# pts panels<sup>®</sup>

## HDL Cholesterol Controls Range Card

		Level 1		Level 2	
Lot #:	HC29	mg/dL	mmol/L	mg/dL	mmol/L
<b>Expiration Date:</b> 2021-12-31					
<b>Single Analyte Test Strips</b>					
HDL Cholesterol* <i>(CardioChek<sup>®</sup> PA analyzers and CardioChek<sup>®</sup> Plus analyzers version 1.08 or below)</i>		≤25-45	≤0.65-1.17	30-65	0.78-1.68
<b>Multi-Analyte Test Strips</b>					
HDL Cholesterol* <i>(CardioChek<sup>®</sup> PA analyzers and CardioChek<sup>®</sup> Plus analyzers version 1.08 or below)</i>		≤15-50	≤0.39-1.30	25-65	0.65-1.68

\* For CardioChek, CardioChek PA, and CardioChek Plus (version 1.08 or below) analyzers.

PTS Diagnostics Customer Service: 877-870-5610 | Fax: 317-870-5608 | Outside of US: +1-317-870-5610  
Email: customerservice@ptsdiagnostics.com