

<b>pts panels<sup>®</sup></b> Multi-Chemistry Controls      Range Card	Level 1		Level 2	
	mg/dL	mmoL/L	mg/dL	mmoL/L
<b>Lot #:</b> <b>MC38</b> <b>Expiration Date:</b> <b>2021-12-31</b>				
<b>Single-Analyte Test Strips</b>				
Cholesterol	110-220	2.85-5.70	140-320	3.63-8.29
HDL Cholesterol* <i>(CardioChek<sup>®</sup> Plus analyzer version 1.09 or above)</i>	≤25-60	≤0.65-1.55	30-80	0.78-2.07
Triglycerides	80-200	0.90-2.26	150-300	1.70-3.39
Glucose/eGLU	40-130	2.22-7.22	200-480	11.10-26.64
<b>Multi-Analyte Test Strips</b>				
Cholesterol	105-225	2.72-5.83	130-320	3.37-8.29
HDL Cholesterol* <i>(CardioChek<sup>®</sup> Plus analyzer version 1.09 or above)</i>	≤15-60	≤0.39-1.55	25-85	0.65-2.20
Triglycerides	70-210	0.79-2.37	140-320	1.58-3.62
Glucose	20-125	1.11-6.94	130-400	7.22-22.20
* For CardioChek, CardioChek PA, and CardioChek Plus (version 1.08 or below) analyzers, HDL Cholesterol requires a separate HDL Cholesterol Control (Cat. No. 722).				
PTS Diagnostics Customer Service: 877-870-5610   Fax: 317-870-5608   Outside of US: +1-317-870-5610 Email: customerservice@ptsdiagnostics.com				
PS-002519 Rev. 11				