

pts panels[®]

HDL Cholesterol Controls

Range Card

Level 1

Level 2

Lot #:

HC34

Expiration Date:

2023-07-31

mg/dL

mmol/L

mg/dL

mmol/L

Single-Analyte Test Strips

HDL Cholesterol*

≤25-50

≤0.65-1.30

30-60

0.78-1.55

Multi-Analyte Test Strips

HDL Cholesterol*

≤15-50

≤0.39-1.30

25-65

0.65-1.68

* For CardioChek Plus (version 1.08 or below), CardioChek PA, and CardioChek analyzers.

PTS Diagnostics Customer Service: 877-870-5610 | Fax: 317-870-5608 | Outside of US: +1-317-870-5610

Email: customerservice@ptsdiagnostics.com

PS-002516 Rev. 9