

| pts panels[®] Multi-Chemistry Controls Range Card | | Level 1 | | Level 2 | |
|---|-------------------|--------------|---------------|--------------|---------------|
| Lot #: | MC42 | mg/dL | mmol/L | mg/dL | mmol/L |
| Expiration Date: | 2023-05-31 | | | | |
| Single-Analyte Test Strips | | | | | |
| Cholesterol | | 130-230 | 3.37-5.96 | 160-310 | 4.14-8.03 |
| HDL Cholesterol* <i>(CardioChek[®] Plus analyzer version 1.09 or above)</i> | | ≤25-60 | ≤0.65-1.55 | 30-70 | 0.78-1.81 |
| Triglycerides | | 90-195 | 1.02-2.20 | 165-300 | 1.86-3.39 |
| Glucose/eGLU | | 40-120 | 2.22-6.66 | 185-400 | 10.27-22.20 |
| Multi-Analyte Test Strips | | | | | |
| Cholesterol | | 100-250 | 2.59-6.48 | 140-340 | 3.63-8.81 |
| HDL Cholesterol* <i>(CardioChek[®] Plus analyzer version 1.09 or above)</i> | | ≤15-60 | ≤0.39-1.55 | 25-80 | 0.65-2.07 |
| Triglycerides | | 70-190 | 0.79-2.15 | 140-300 | 1.58-3.39 |
| Glucose | | 25-130 | 1.39-7.22 | 120-400 | 6.66-22.20 |
| *HDL Cholesterol requires separate HDL Cholesterol Controls (REF 722) for CardioChek Plus (version 1.08 or below), CardioChek PA, and CardioChek analyzers. | | | | | |
| PTS Diagnostics Customer Service: 877-870-5610 Fax: 317-870-5608 Outside of US: +1-317-870-5610 Email: customerservice@ptsdiagnostics.com | | | | | |