

pts panels[®] HDL Cholesterol Controls Range Card		Level 1		Level 2	
Lot #:	HC35	mg/dL	mmol/L	mg/dL	mmol/L
Expiration Date:	2023-09-30				
Single-Analyte Test Strips					
HDL Cholesterol*		≤25-45	≤0.65-1.17	30-60	0.78-1.55
Multi-Analyte Test Strips					
HDL Cholesterol*		≤15-50	≤0.39-1.30	30-65	0.78-1.68
<p>* For CardioChek Plus (version 1.08 or below), CardioChek PA, and CardioChek analyzers.</p> <p>PTS Diagnostics Customer Service: 877-870-5610 Fax: 317-870-5608 Outside of US: +1-317-870-5610 Email: customerservice@ptsdiagnostics.com</p>					
					PS-002516 Rev. 9