

pts panels[®] Multi-Chemistry Controls Range Card		Level 1		Level 2	
		mg/dL	mmol/L	mg/dL	mmol/L
Lot #:	MC43				
Expiration Date:	2023-10-31				
Single-Analyte Test Strips					
Cholesterol		130-210	3.37-5.44	150-290	3.89-7.51
HDL Cholesterol* <i>(CardioChek[®] Plus analyzer version 1.09 or above)</i>		≤25-55	≤0.65-1.42	30-70	0.78-1.81
Triglycerides		100-185	1.13-2.09	160-280	1.81-3.16
Glucose/eGLU		40-120	2.22-6.66	195-350	10.82-19.43
Multi-Analyte Test Strips					
Cholesterol		100-230	2.59-5.96	150-320	3.89-8.29
HDL Cholesterol* <i>(CardioChek[®] Plus analyzer version 1.09 or above)</i>		≤15-55	≤0.39-1.42	25-80	0.65-2.07
Triglycerides		80-190	0.90-2.15	140-300	1.58-3.39
Glucose		30-130	1.67-7.22	150-380	8.33-21.09
* HDL Cholesterol requires separate HDL Cholesterol Controls (REF 722) for CardioChek Plus (version 1.08 or below), CardioChek PA, and CardioChek analyzers.					
PTS Diagnostics Customer Service: 877-870-5610 Fax: 317-870-5608 Outside of US: +1-317-870-5610 Email: customerservice@ptsdiagnostics.com					
PS-002519 Rev. 15					