

pts panels[®] Multi-Chemistry Controls Range Card	Level 1		Level 2	
	mg/dL	mmol/L	mg/dL	mmol/L
Lot #: MC44 Expiration Date: 2023-10-31				
Single-Analyte Test Strips				
Cholesterol	120-200	3.11-5.18	170-310	4.40-8.03
HDL Cholesterol* <i>(CardioChek[®] Plus analyzer version 1.09 or above)</i>	≤25-55	≤0.65-1.42	30-80	0.78-2.07
Triglycerides	100-180	1.13-2.03	170-290	1.92-3.28
Glucose/eGLU	40-110	2.22-6.11	200-340	11.10-18.87
Multi-Analyte Test Strips				
Cholesterol	110-230	2.85-5.96	150-345	3.89-8.94
HDL Cholesterol* <i>(CardioChek[®] Plus analyzer version 1.09 or above)</i>	≤15-65	≤0.39-1.68	25-90	0.65-2.33
Triglycerides	90-200	1.02-2.26	150-300	1.70-3.39
Glucose	40-120	2.22-6.66	170-370	9.44-20.54
* HDL Cholesterol requires separate HDL Cholesterol Controls (REF 722) for CardioChek Plus (version 1.08 or below), CardioChek PA, and CardioChek analyzers.				
PTS Diagnostics Customer Service: 877-870-5610 Fax: 317-870-5608 Outside of US: +1-317-870-5610 Email: customerservice@ptsdiagnostics.com				
				PS-002519 Rev. 15