

# pts panels<sup>®</sup>

Multi-Chemistry Controls Range Card

Level 1

Level 2

Lot #: MC45  
Expiration Date: 2023-11-30

mg/dL

mmol/L

mg/dL

mmol/L

## Single-Analyte Test Strips

Cholesterol	130-210	3.37-5.44	170-300	4.40-7.77
HDL Cholesterol* <i>(CardioChek<sup>®</sup> Plus analyzer version 1.09 or above)</i>	≤25-55	≤0.65-1.42	30-80	0.78-2.07
Triglycerides	100-180	1.13-2.03	170-290	1.92-3.28
Glucose/eGLU	40-110	2.22-6.11	200-340	11.10-18.87

## Multi-Analyte Test Strips

Cholesterol	110-230	2.85-5.96	150-350	3.89-9.07
HDL Cholesterol* <i>(CardioChek<sup>®</sup> Plus analyzer version 1.09 or above)</i>	≤15-60	≤0.39-1.55	25-90	0.65-2.33
Triglycerides	90-180	1.02-2.03	160-300	1.81-3.39
Glucose	40-120	2.22-6.66	170-360	9.44-19.98

\* HDL Cholesterol requires separate HDL Cholesterol Controls (REF 722) for CardioChek Plus (version 1.08 or below), CardioChek PA, and CardioChek analyzers.

PTS Diagnostics Customer Service: 877-870-5610 | Fax: 317-870-5608 | Outside of US: +1-317-870-5610  
Email: customerservice@ptsdiagnostics.com

PS-002519 Rev. 15