

<b>pts panels<sup>®</sup></b> Multi-Chemistry Controls    Range Card		Level 1		Level 2	
<b>Lot #:</b>	<b>MC46</b>	<b>mg/dL</b>	<b>mmol/L</b>	<b>mg/dL</b>	<b>mmol/L</b>
<b>Expiration Date:</b>	<b>2024-03-31</b>				
<b>Single-Analyte Test Strips</b>					
Cholesterol		120-200	3.11-5.18	170-280	4.40-7.25
HDL Cholesterol* <i>(CardioChek<sup>®</sup> Plus analyzer version 1.09 or above)</i>		≤25-55	≤0.65-1.42	30-65	0.78-1.68
Triglycerides		110-180	1.24-2.03	160-270	1.81-3.05
Glucose/eGLU		40-110	2.22-6.11	200-350	11.10-19.43
<b>Multi-Analyte Test Strips</b>					
Cholesterol		100-220	2.59-5.70	140-350	3.63-9.07
HDL Cholesterol* <i>(CardioChek<sup>®</sup> Plus analyzer version 1.09 or above)</i>		≤15-55	≤0.39-1.42	25-80	0.65-2.07
Triglycerides		80-170	0.90-1.92	160-270	1.81-3.05
Glucose		40-115	2.22-6.38	130-350	7.22-19.43
* HDL Cholesterol requires separate HDL Cholesterol Controls (REF 722) for CardioChek Plus (version 1.08 or below), CardioChek PA, and CardioChek analyzers.					
PTS Diagnostics Customer Service: 877-870-5610   Fax: 317-870-5608   Outside of US: +1-317-870-5610 Email: customerservice@ptsdiagnostics.com					
PS-002519 Rev. 15					