

**pts panels<sup>®</sup>**

HDL Cholesterol Controls

Range Card

Level 1

Level 2

Lot #:

HC37

Expiration Date:

2025-05-31

mg/dL

mmol/L

mg/dL

mmol/L

**Single-Analyte Test Strips**

HDL Cholesterol\*

≤25-45

≤0.65-1.17

30-60

0.78-1.55

**Multi-Analyte Test Strips**

HDL Cholesterol\*

≤15-50

≤0.39-1.30

25-65

0.65-1.68

\* For CardioChek Plus (version 1.08 or below), CardioChek PA, and CardioChek analyzers.

PTS Diagnostics Customer Service: 877-870-5610 | Fax: 317-870-5608 | Outside of US: +1-317-870-5610

Email: [customerservice@ptsdiagnostics.com](mailto:customerservice@ptsdiagnostics.com)

PS-002516 Rev. 9