

pts panels[®]

Multi-Chemistry Controls

Range Card

Level 1

Level 2

Lot #: MC50

Expiration Date: 2025-07-31

mg/dL

mmol/L

mg/dL

mmol/L

Single-Analyte Test Strips

Cholesterol	120-210	3.11-5.44	150-310	3.89-8.03
HDL Cholesterol* <i>(CardioChek[®] Plus analyzer version 1.09 or above)</i>	≤25-50	≤0.65-1.30	30-65	0.78-1.68
Triglycerides	90-180	1.02-2.03	160-260	1.81-2.94
Glucose/eGLU	30-110	1.67-6.11	160-360	8.88-19.98

Multi-Analyte Test Strips

Cholesterol	100-220	2.59-5.70	130-320	3.37-8.29
HDL Cholesterol* <i>(CardioChek[®] Plus analyzer version 1.09 or above)</i>	≤15-50	≤0.39-1.30	25-75	0.65-1.94
Triglycerides	70-170	0.79-1.92	140-280	1.58-3.16
Glucose	40-115	2.22-6.38	120-340	6.66-18.87

* HDL Cholesterol requires separate HDL Cholesterol Controls (REF 722) for CardioChek Plus (version 1.08 or below), CardioChek PA, and CardioChek analyzers.

PTS Diagnostics Customer Service: 877-870-5610 | Fax: 317-870-5608 | Outside of US: +1-317-870-5610

Email: customerservice@ptsdiagnostics.com

PS-002519 Rev. 15